ExCell

EIV22024E

ExCell Bio

ResiQuant[®] CHO HCP ELISA Kit-FAST

For Research and Manufacturing Use

Not Intended for Diagnostic and Therapeutic Use

User Manual

Catalog Number

CRH00-3051S

CRH00-3051

CRH00-3052




PRODUCT DESCRIPTION

The expression of therapeutic proteins in the CHO cells lines is economical and convenient for commercial drug production. However, the host cell proteins (HCPs) from CHO cells as a type of concomitant impurities are one of the major concerns in terms of drug safety. These impurities may diminish the efficacy of drug therapy and cause side effects or unexpected immune reactions. Therefore, control and monitor of the residual CHO HCP impurities is world widely required by all corresponding regulatory authorities. This product is a double-antibody sandwich ELISA kit. The antibodies were generated from rabbits and goats against the antigens harvested from the fermentation supernatant of CHO cells, cultured in the medium free of proteins. This kit can be used throughout the manufacturing process and for product release. For the first time, it is recommended to conduct a suitability study to confirm the sample matrix effect and determine the appropriate dilution conditions for sample testing.

PRINCIPLES OF THE TEST

The assay employs a sandwich immunoassay methodology for the quantitation of CHO HCP. Microplates are coated with antibodies that exhibit specificity for CHO HCP. Anti-CHO HCP antibodies labeled by HRP and standards or samples are added to the wells simultaneously, during incubation, CHO HCP in standard and samples forms immune complexes with antibodies pre-coated and labeled. After washing away unbound substances, the chromogenic substrate 3,3',5,5'-tetramethylbenzidine (TMB) is introduced, initiating an enzymatic reaction. This reaction involves HRP-mediated oxidation of TMB in the presence of hydrogen peroxide (H₂O₂), yielding a blue-colored product that exhibits maximum absorbance at 655 nm. The enzymatic reaction is subsequently halted by the addition of a stop solution, which turns the color to yellow with peak absorbance at 450 nm. The optical density (OD) value at 450 nm, measured by an enzyme-linked immunosorbent assay (ELISA) reader, is proportional to the concentration of CHO HCP within the samples. Quantitative determination of CHO HCP is achieved with a simultaneously prepared standard curve, using the standard supplied within the kit.

Schematic diagram :



Web: www.excellbio.com

| PRODUCT PERFORMANCE

ExCell

1. Sensitivity: The limit of detection (LOD) for CHO HCP is 0.5 ng/mL, and the low limit of quantification (LLOQ) is 1.0 ng/mL.

2. Precision: The coefficient of variation (CV) for sample concentrations within the quantitation range is no more than 20% for both intra and inter-assay.

3. Specificity: No cross-reactivity has been observed with *E.coli*, Vero, HEK293 whole cell lysates.

PRODUCT APPLICATION

This product is a general assay kit designed for the quantitative detection of CHO HCP impurities in samples.

| PRODUCT SPECIFICATIONS

| Catalog Number | Name | Size |
|----------------|--|------|
| CRH00-3051S | ResiQuant [®] CHO HCP ELISA Kit -FAST | 48T |
| CRH00-3051 | ResiQuant [®] CHO HCP ELISA Kit -FAST | 48T |
| CRH00-3052 | ResiQuant [®] CHO HCP ELISA Kit -FAST | 96T |

REAGENTS PROVIDED

| Name | 96 Tests | 48 Tests | Store | | |
|---|----------|----------|---------------------------------|--|--|
| CHO HCP 2G Microplate | 8×12 | 8× 6 | 2°C to 8°C | | |
| CHO HCP Standard Standards at 250, 100, 40, 16, 6.4, 2.6, 1 and 0 ng/mL, 1 mL/vial | 8 | 8 | 2°C to 8°C | | |
| 100× CHO HCP HRP-Antibody | 60 µL | 30 µL | 2°C to 8°C | | |
| CHO HCP Assay Diluent | 2× 25 mL | 25 mL | 2°C to 8°C | | |
| 20× Wash Buffer Concentrate | 30 mL | 30 mL | 2°C to 8°C | | |
| Substrate Solution | 12 mL | 6 mL | 2°C to 8°C (Light-Sensitive) | | |
| Stop Solution | 12 mL | 12 mL | 2°C to 8°C | | |
| Plate Sealer | 3 | 2 | / | | |

Tel: (+86) 4008205021

INSTRUCTION FOR USE

I. Materials required but not provided

1. Microplate reader capable of measuring absorbance at 450 nm (If wavelength correction is needed, additional wavelength at 630 nm or 570 nm is required).

1. 0.5-10 µL, 2-20 µL, 20-200 µL, 100-1000 µL adjustable micropipettes with disposable tips.

- 2. Microtiter plate rotator.
- 3. Deionized water.

ExCell

II. Sample collection

1. Samples should be clarified and any precipitate should be removed by centrifugation.

2. The sample should be diluted in appropriate proportions according to the condition of the sample. (It is recommended to complete the applicability study of the kit to confirm the appropriate dilution conditions for sample detection).

III. Preparation of reagents

1. It is recommended to remove the assay kit from the refrigerator 20 minutes in advance to allow it to equilibrate to room temperature.

1. Dilute $20 \times$ Wash Buffer Concentrate with deionized or distilled water to prepare the wash buffer. Store the unused solution at 2°C to 8°C.

2. CHO HCP HRP-Antibody Working Solution: Freshly prepare every time the working solution by diluting $100 \times$ CHO HCP HRP-Antibody with CHO HCP Assay Diluent to $1 \times$ in an amount sufficient for the current experiment.

IV. Washing

1. Automatic washing machine: set the volume to 300 μ L, and the interval between aspiration and dispensing is 10 seconds. Wash the plate 5 times.

2. Manual plate washing: add 300 μ L of wash buffer into each well, allow it to sit for 10 seconds, then shake off the liquid in the well. Invert the plate and blot it against clean paper towels. Wash the plate 5 times.

V. Assay procedure

1. Bring all reagents and samples to room temperature before use. Remove excess microplate strips from the plate frame, return them to the foil pouch containing a desiccant pack, reseal tightly and store at 2°C to 8°C.

2. Set the blank well (if the plate is measured by dual-wavelength, the blank may not be necessary).

3. Prepare samples and standards in advance.

4. Add 100 μ L per well of either standards or samples, then add 50 μ L per well of CHO HCP HRP-Antibody Working Solution immediately. Seal the wells with the plate sealer. Incubate at 25°C and shake for 60 minutes, using a microplate shaker (500 rpm).

- 5. Remove the liquid from the wells and wash the plate 5 times.
- 6. Add 100 μL of Substrate Solution to the wells, Incubate at 20°C to 25°C in the dark for 15 minutes.
- 7. Add 100 μ L of Stop Solution to well.

8. Read the optical density (OD) of each well within 10 minutes using a microplate reader set to 450 nm.


If wavelength correction is available, set to 570 nm or 630 nm. For single wavelength reading, the OD value of

Web: www.excellbio.com Tel: (+86) 4008205021 Email: marketing@excellbio.com



blank well shall be subtracted from that of each Standard or sample.

VI. Assay Procedure Summary



ANALYSES OF RESULTS

ExCell

1. Create a standard curve by plotting the mean optical density (OD) of each standard concentration as the y-axis against the standard concentration on the x-axis. It is recommended to use a four-parameter fitting equation for the standard curve. Calculate the concentration of samples based on their OD values.

2. The coefficient of determination (\mathbb{R}^2) for the standard curve should be no less than 0.99. Discard those obvious abnormal data, the deviation between the back calculated and the theoretical concentration of each standard should be within ±20%, (±25% for the upper or lower limits of quantitation).

3. The concentration of samples should be calculated with the standard curve in the parallel assay.

4. If the OD value of the sample exceeds the high end of the standard curve, the sample should be appropriately diluted before further measurement. Be aware that the dilution factor should be counted with back calculated concentration together.

| Standard | OD 450-630 | | Concentration (ng/mL) | | | 3 | | | | | ٦ | |
|-----------------------|-------------------|-------|-----------------------|-------|-------|---------|---------------------|----------|----------|----------------------|-------|-----|
| (ng/mL) | 1 | 2 | Average | 1 | 2 | Average | 2.5 - | | | ļ | | |
| 250 | 2.606 | 2.493 | 2.550 | 261.0 | 239.2 | 250.1 | | | | | | |
| 100 | 1.515 | 1.465 | 1.490 | 102.9 | 97.8 | 100.3 | 2 - | | | | | |
| 40 | 0.780 | 0.770 | 0.775 | 40.0 | 39.3 | 39.6 | 20 | | | | | |
| 16 | 0.407 | 0.399 | 0.403 | 16.3 | 15.9 | 16.1 | 051 05 00 | | / | | | |
| 6.4 | 0.217 | 0.231 | 0.224 | 6.2 | 6.9 | 6.6 | 1- | | | | Ļ | |
| 2.6 | 0.140 | 0.143 | 0.142 | 2.6 | 2.7 | 2.6 | | ø | | | | |
| 1.0 | 0.104 | 0.104 | 0.104 | 1.0 | 1.0 | 1.0 | 0.5 - | | | | | |
| 0 | 0.076 | 0.074 | 0.075 | | | | | 7 | | | | |
| R ² | 0.99999 | | | | | | | , | 50 50 | 100 entration (no | 200 | 250 |

Representative Standard Curve:

Note: This graph is provided for reference purpose only. Data analysis was using ELISA Calc software with a four-parameter fitting model.

1. Store the reagents at 2°C to 8°C. The reconstituted standard stock solution should be kept at 18°C or below.

2. The concentrated Antibody-HRP is supplied in a small volume. Turbulence or inversion happened during transportation may cause spray of the solution elsewhere on the wall or the cap. Therefore, please spin the tubes shortly before use to harvest all liquid to the bottom of the tube.

3. In case crystal sediment appears in the $20 \times$ Wash Buffer Concentrate, warm it up at 37°C to dissolve the sediment before use.

- 4. Do not mix or substitute reagents with those from other lots or different kits.
- 5. Ensure thoroughly mixing in preparing the solutions to ensure the reactions are consistent.
- 6. It is recommended to test all samples and standards in duplicate form at least.

DISCLAIMER

ExCell

NOTES

1. The product should be used according to the instructions in the manual. If the experimenter fails to operate according to the instructions, our company will not be responsible for any deviation in product performance caused by this.

2. The product is only used for scientific research and commercial production, and is not suitable for clinical diagnosis and treatment. Otherwise, all consequences arising shall be borne by the experimenter, and our company shall not be responsible.